



SNAIGĖ MAX line PHARMACEUTICAL REFRIGERATORS ENSURING OPTIMAL MEDICATION STORAGE

Certified DIN 13277

MAX line refrigerators have all the essential features and performance functions required to meet the performance criteria of the DIN 13277 standard. This standard provides the requirements and test specifications for medical and laboratory refrigerators and freezers.

Test method EN 60068-3

SNAIGĖ medical refrigerators are tested according to EN 60068-3 in order to ensure the stability and constancy of the internal temperature. The test method involves measuring the temperature at 15 different points inside the refrigerator.

Ecological and environmentally friendly

MAX line refrigerators are manufactured using only non-harmful, non-ozone-depleting and non-greenhouse-effect refrigerating and insulating materials. Natural resource-saving materials and electricity from renewable energy sources are used in the refrigerator manufacturing process.

MD40DD-P301SEQ / MD40DD-P600SEQ

MAX line Pharmacy refrigerator

KEY FEATURES



Electronic control system

An electronic controller installed in the MAX line refrigerators provides continuous monitoring of the temperature inside the refrigerator.



Temperature stability

The refrigerator is designed to maintain maximum stable temperature and minimal fluctuations. The maximum possible temperature change within the entire internal space of the refrigerator is only ± 2.9 °C. Thus, the internal microclimate is maintained, creating suitable conditions for the storage of medications and medicines.



Cooling system with circulating air

The SNAIGĖ smart fan for medical refrigerators operates very efficiently. It turns off as soon as the refrigerator door is opened, so the air circulation in the refrigerator stops. Once the door is closed, the fan starts again and the air circulation is resumed. This fan operation prevents rapid temperature changes and helps to maintain a stable temperature inside the refrigerator, contributing to the preservation of medicines and medications.



Protection against contact with cold surfaces

The evaporator of the refrigerator has a special cover so that medications and medicines do not have direct contact with the evaporator and do not get cold. In addition, this protection against contact with cold surfaces increases the temperature stability.



MD40DD-P600SEQ

MD40DD-P301SEQ

KEY FEATURES



Integrated data logger

The memory of the MD40DD-P301SEQ electronic controller stores the last 100 alarm records with the date and duration of the alarm, as well as the temperature measurement summary with the measurement date and time. These saved data can be transferred to a USB storage device.



Warning of electricity supply failure

In the event of the electricity supply failing, the battery of the MD40DD-P301SEQ provides sufficient energy for the registration of temperature readings and the operation of the warning signal (audio and visual) for the next 12 hours.



Warning temperature alarm

In order to ensure the maximum safety of stored medicines and medications, the integrated warning temperature alarm sends a warning signal with sound and image when a certain temperature value exceeds or falls below the set limit.



Warning register

Warning signals must be monitored to ensure the stored medicines or medications retain their quality and integrity. Important data are recorded in the register: date of warning, duration, highest temperature. The memory of the MD40DD-P301SEQ electronic controller stores the last 100 warning records.



+2 °C temperature limit

If the temperature drops below +2 °C for some critical reasons, the electronic controller immediately activates additional temperature maintenance functions along with a warning signal in order to protect the temperature-sensitive medicines and medications stored in the refrigerator.



LAN interface

SNAIGÉ medical refrigerators can be connected to a LAN network (using external systems) through the integrated LAN interface RS485. It is convenient to monitor, control and digitally document the temperature of medicines stored in the refrigerator.



Access to an additional external sensor

An additional external sensor can be attached if required for additional control of the temperature of medication storage in the refrigerator. SNAIGÉ medical refrigerators are equipped with a 12.5 mm diameter access for an additional external sensor, which helps avoid any negative effect of the cables of the external sensor on the sealing strip of the door, the tightness of the refrigerator and the risk to the personnel caused by the cables of the sensor.



Self-closing door

The door of the SNAIGÉ medical refrigerator closes automatically, without any risk of the door bouncing. This ensures convenience and safety for medical personnel in their daily use of the refrigerator, and also completely eliminates the risk of the refrigerator door being accidentally left open. The self-closing door contributes to maintaining a stable temperature, as well as supporting the integrity of stored medicines and medications.



TECHNICAL DATA

Pharmacy refrigerator	MD40DD-P301SEQ MXXXXX1SN4B	MD40DD-P600SEQ MXXXXX1SN4B
EAN code	4770104040092	4770104040054
Gross capacity, l	385	
Net capacity, l	290	
Factory temperature setting, °C	+5	
Fluctuation, °C	0,5	
Gradient, °C	2,9	
Increase in temperature from +5 °C to +10 °C (empty, +25 °C TU)	40 min.	
Decrease in temperature from +25 °C to +5 °C (empty, +25 °C TU)	40 min.	
Ambient temperature range, °C	from +10 to +35	
Airborne acoustical noise emissions, dB(A) re1 pW	50	
Refrigerant	R600a	
Daily consumption, kWh/24h	1,903	1,440
Annual energy consumption, kWh/year	695	525
Rated voltage, V	220 - 240	
Rated frequency, Hz	50	
Access port	Yes	
Cooling system	Fan forced cooling	
Defrost	Automatic	
Volt-free alarm contact	Yes	
Interior light	Ceiling LED light	
Number of shelves	5 adjustable +1 bottom	
Shelf size (W x D), mm	380 x 495, bottom 445 x 180	
Maximum load per shelf, kg	40	
Shelf material	Wire shelves coated by white polymer paint	
Side wall material	Steel, white painted	
Inner cabinet material	HIPS polystyrene	
Type of door	Double panel selective glass in aluminum profile	Full-panel door, white painted
Reversible door	Yes	
Self closing door	Yes	
Type of lock	Mechanical	
Adjustable foot at the front	2	
Wheels behind	2	
External dimensions (H x W x D), mm	2025 x 600 x 620	
Internal dimensions (H x W x D), mm	1645 x 500 x 418	
Package dimensions (H x W x D), mm	2080 x 620 x 630	
Weight net, kg	80,5	
Weight gross, kg	86	

Contact us: snaige@snaige.lt

- Phone +370 315 56206
- Mob.: +370 656 35923
- AB SNAIGĖ, Pramonės str.
6 LT-62175, Alytus,
Lithuania

More information:

www.snaigeprofessional.lt

